

pH Meter 3 point

Microprocessor based

Model: CGPH-100



SALIENT FEATURES:

- Large back light LCD display shows mode indicator pH value, temperature
 - Up to 3 point push -button calibration with automatic buffer recognition
 - Selectable pH buffer standards and temperature units.
 - Automatically calculates and display electrode slope after calibration
 - Automatic temperature compensation (ATC) ensure highly accurate measurements over the entire range.
 - Automatic endpoint function freezes the stable measured value
- Setup menu allows setting the pH buffer group, number of calibration points, temperature unit, auto-power off, etc.
 - Reset function can resumes all setting to factory default value

Galaxy PH-100 is an easy-to-use benchtop pH meter designed for daily testing and quality control.

This meter has a clear and bright back light LCD display that show the measured pH value, temperature, mode indicators.

Digital pH Meter

